



General order information

Material

- StandardX-light
- Lead free

All three materials are in generally very flexible and robust. The difference is only in weight. The X-light material is approximately 20% lighter and the lead-free one is roughly 25% lighter compared with the standard material.

The protection level is for all materials the same and belongs to the lead equivalent information of 0.25 / 0.35 / 0.50 Pb or more. Wiroma aprons are CE tested for personal protection and are checked between 50 an 150 kV. Of course, all our products can be used as well for patient protection.

Our test centre for CE certificate is TÜV SÜD Product Service GmbH in Munich with the ID number 0123. The requirements conform to article 11, B, No. 1(c) of Council Directive 89/686/EEC for personal protective equipment. Our products are already tested to new Norm IEC 61331-3:2002 and 61331-1:2014. Our products have also passed the CE-Type Approval.

Sizes


Aprons are available in different sizes; see the price list. Please consult our measurements table below so that we can deliver protective wear in the correct size.

FRONTSCHÜRZE TABLIER FRONTAL FRONT APRON

Breite = von Schulter zu Schulter (gemäss Abbildung)

Largeur = d'épaule à épaule (comme image)

Width = from shoulder to shoulder (as drawing)



-ongueur : de la moitié des épaules à la moitié des genoux

-ength: measure mid-shoulder to mid-knee

MANTELSCHÜRZE TABLIER-MANTEAU WRAP AROUND APRON

Brust- oder Bauchumfang. Nur das grössere Mass angeben.

Tour de la poitrine ou du ventre. A retenir que la mesure la plus grande.


Measure circumference of chest and abdomen. State only the larger of the two measurements.

Langueur : de la moitié des épaules à la moitié des genoux Length: measure mid-shoulder to mid-knee. TOP UND JUPE VESTE ET JUPE TOP AND SKIRT

Top: Brust- oder Bauchumfang. Nur das grössere Mass angeben. Jupe: Taille

Top: la circonférence de la poitrine ou du ventre. A retenir que la mesure la plus grande. Jupe: taille

Top: measure circumfrence of chest and abdomen. State only the larger of the two measurements. Skirt: measure waist



Farbcode-Nr. / Nr. De code du coloirs / Colour codes

Hauptfarbe und Farbe der Brusttasche / Couleur principale et couleur de la poche / Main colour and colour of the brest pocket

Material Nylon/Matière Nylon/Material nylon



petrol dunkelblau/bleu foncé/dark blue blau/bleu/blue dunkelrot/rouge foncé/dark red türkis/turquoise/turquoise

Material Magic/Matière Magic/Material Magic



"Magic"

- Flexible, elasticated and very hard-wearing material. Betadine, Hexaquart plus, blood, contrast media and, to some extent, aggressive fluids can be easily removed with warm soapy water.
- Can be used in surgery, angiography and in all work areas that have an increased risk of contamination.

"Nylon"

- Mechanically-speaking, it is more robust and more scratch-resistant than Magic.
- Can be used in general medicine and in veterinary practice.

Monograms

- Up to 21 letters per line possible
- Example: 'Prof. Hansrudolf Müller' 'Radiology'

WIROMA Options

- Inner chest pocket for dosimeter
- Affixed logos and monograms
- Transparent insert pocket for name card
- Additional sleeves
- Thyroid collar either sown on or with press studs
- Children's motifs
- Additional fixing ties (Velcro)
- All radiation protection aprons can also be made in special sizes

WIROMA USP

- Integrated chest pocket (removable)
- Resistant to commercial disinfectants
 (Hexaquart plus, etc.)
- RFID chip controlled

WIROMA by Lanz-Anliker AG	
	Clothing
Your specialist for radiation protection	Front Aprons

Length XSV 90 XSS 100 XSM 110 SXL 120 XSE 130 SV 90 SS 100 SM 110

SL 120 SE 130 MV 90 MS 100 MM 110

ML 120 ME 130 LV 90 LS 100 LM 110

LL 120 LE 130

413-AWC	
WIROMA Pb 0.25 mm	X-SMALL Width 50 cm < 99 cm
	SMALL Width 60 cm Up to ca. 115 cm
	MEDIUM Width 65 cm Up to ca. 125 cm
	LARGE Width 75 cm Up to ca. 140 cm

Optimal weight distribution is achieved with the AWC apron thanks to the special form.

The front apron is also available with an integrated thyroid collar.

Option: Additional velcro tape to fix the shoulder bands.



	Length
	XSV 90
X-SMALL	XSS 100
Width 50 cm	XSM 110
< 99 cm	SXL 120
	XSE 130
	SV 90
SMALL	SS 100
Width 60 cm	SM 110
Up to ca. 115 cm	SL 120
	SE 130
	MV 90
MEDIUM	MS 100
Width 65 cm	MM 110
Up to ca. 125 cm	ML 120
	ME 130
	LV 90
LARGE	LS 100
Width 75 cm	LM 110
Up to ca. 140 cm	LL 120
	LE 130

The WBT apron provides great comfort and is very adaptable due to a wide cross-over belt that has a practical plastic buckle. Less contamination (veterinary surgeon) and very hygienic.

The front apron is also available with an integrated thyroid collar.

Option: Additional velcro tape to fix the shoulder bands.

Operating aprons

421-SDO



	Length
	XSV 90
X-SMALL	XSS 100
Width 50 cm	XSM 110
< 99 cm	SXL 120
	XSE 130
	SV 90
SMALL	SS 100
Width 60 cm	SM 110
Up to ca. 115 cm	SL 120
	SE 130
	MV 90
MEDIUM	MS 100
Width 65 cm	MM 110
Up to ca. 125 cm	ML 120
	ME 130
	LV 90
LARGE	LS 100
Width 75 cm	LM 110
Up to ca. 140 cm	LL 120
	LE 130

As soon as there is no more radiation in the operating theatre, these aprons can be removed from underneath the operating clothing by undoing the Velcro straps.

For the rest of the time in the operating theatre you have no need to wear any unnecessary weight.

WIROMA by Lanz-Anliker AG	
	Clothing
Your specialist for radiation protection	Coat Aprons

431-LSW



	Length
	XSV 90
X-SMALL < 99 cm	XSS 100
	XSM 110
	SXL 120
	XSE 130
	SV 90
CMALL	SS 100
SMALL ~ 100 cm	SM 110
~ 100 cm	SL 120
	SE 130
	MV 90
MEDIUM	MS 100
	MM 110
~ 110 cm	ML 120
	ME 130
	LV 90
LARGE	LS 100
~ 120 cm	LM 110
~ 120 Cm	LL 120
	LE 130
	XLV 90
X-LARGE	XLS 100
~ 130 cm	XLM 110
~ 130 cm	XLL 120
	XLE 130

Thanks to its wide stretch belt, the 431-LSW coat apron takes the weight off the shoulders. This model provides optimum shoulder movement and the highest comfort.

The height of the exterior belt is adjustable.

431-LSW GI (with inside belt)



Length XSV 90 XSS 100 XSS 100 XSM 110 99 cm SXL 120 XSE 130 SV 90 SMALL ~ 100 cm SE 130 SK 110 SS 100 SMALL ~ 100 cm MEDIUM MEDIUM MS 100 MM 110 NL 120 LARGE LV 90 LARGE ~ 120 cm LL 120 LE 130 XLV 90 XLS 100 XLL 120 XLL 120 XLE 130		
X-SMALL < 99 cm XSS 100 XSM 110 SXL 120 XSE 130 SV 90 SSMALL ~ 100 cm SS 100 SM 110 SL 120 SE 130 MV 90 MS 100 MM 110 ML 120 ME 130 LV 90 LARGE ~ 120 cm LARGE ~ 120 cm LX 90 LX 120 MK 110 ML 120 ME 130 LV 90 LX 120 LX 120 ME 130 LV 90 LX 120 LX 120 XLX 12		Length
X-SMALL XSM 110 < 99 cm XSM 110 SXL 120 XSE 130 SV 90 SMALL SM 110 SV 90 SS 100 SS 100 SM 110 SL 120 SE 130 MV 90 MS 100 MM 110 ML 120 ME 130 LV 90 LS 100 LV 90 LS 100 LL 120 LS 100 LL 120 LS 100 XLV 90 XLS 100 XLV 90 XLS 100 XLS 100 XLL 120 XLS 100 XLL 120 XLL 120 XL 120	Y-SMALL	XSV 90
 < 99 cm XSM 110 SXL 120 XSE 130 SV 90 SS 100 SS 100 SS 100 SS 100 SE 130 MV 90 MS 100 MM 110 ML 120 MK 100 MM 110 ML 120 LV 90 LS 100 LV 90 LV 90 LS 100 LV 90 LS 100 LV 90 LV 90 LV 90 XLV 90 XLU 120 		XSS 100
SXL 120 XSE 130 SV 90 SS 100 SMALL ~ 100 cm SE 130 SK 120 SE 130 MV 90 MS 100 MK 110 ~ 110 cm ME 130 LV 90 LARGE ~ 120 cm LV 90 LL 120 LE 130 XLV 90 XLS 100 XLS 100 XLM 110 XLL 120	-	XSM 110
SV 90 SMALL SS 100 ~ 100 cm SM 110 SL 120 SE 130 MV 90 MS 100 MM 110 ML 120 ME 130 LV 90 LARGE LV 90 - 120 cm LL 120 X-LARGE XLV 90 X-LARGE XLS 100 - 130 cm XLL 120	< 99 011	SXL 120
SMALL SS 100 ~ 100 cm SM 110 ~ 100 cm SL 120 SE 130 MV 90 MEDIUM MS 100 ~ 110 cm ML 120 ME 130 LV 90 LARGE LM 110 ~ 120 cm LS 100 LL 120 LE 130 XLV 90 XLS 100 XLS 100 XLM 110 XLS 100 XLL 120		XSE 130
SMALL SM 110 ~ 100 cm SL 120 SE 130 MV 90 MEDIUM MV 90 ~ 110 cm MM 110 ML 120 ME 130 LARGE LV 90 LARGE LM 110 ~ 120 cm LL 120 KL 120 KL 100 X-LARGE XLV 90 X-LARGE XLM 110 ~ 130 cm XLL 120		SV 90
~ 100 cm ~ 100 cm SL 120 SE 130 MV 90 MS 100 MM 110 ~ 110 cm LARGE ~ 120 cm X-LARGE ~ 130 cm SL 120 MV 90 LS 100 LV 90 LS 100 LM 110 LL 120 LE 130 XLV 90 XLS 100 XLS 100 XLM 110 XLL 120	CMALL	SS 100
SL 120 SE 130 MV 90 MS 100 MM 110 ~ 110 cm ML 120 ME 130 LV 90 LARGE ~ 120 cm LL 120 LL 120 LL 120 LL 120 LS 100 XLV 90 XLV 90 XLS 100 XLS 100 XLM 110 XLL 120	-	SM 110
MV 90 MEDIUM MS 100 ~ 110 cm MM 110 ML 120 ME 130 LARGE LV 90 ~ 120 cm LM 110 LL 120 LS 100 LL 120 LS 100 X-LARGE XLV 90 X-LARGE XLS 100 XLL 120 XLL 120		SL 120
MEDIUM MS 100 ~ 110 cm MM 110 ML 120 ME 130 LV 90 LS 100 LARGE LM 110 ~ 120 cm LL 120 LE 130 XLV 90 X-LARGE XLS 100 ~ 130 cm XLL 120		SE 130
MEDIUM MM 110 ~ 110 cm ML 120 ME 130 LV 90 LARGE LN 100 ~ 120 cm LL 120 LE 130 XLV 90 X-LARGE XLS 100 ~ 130 cm XLL 120		MV 90
~ 110 cm MM 110 ML 120 ME 130 LV 90 LS 100 LM 110 LM 110 LL 120 LE 130 XLV 90 XLS 100 XLS 100 XLV 90 XLS 100 XLS 100 XLV 110 XLL 120		MS 100
ML 120 ME 130 LV 90 LS 100 LM 110 LL 120 LE 130 X-LARGE ~ 130 cm XLX 90 XLL 120 XLX 100 XLL 120	-	MM 110
LARGE ~ 120 cm X-LARGE ~ 130 cm X-LARGE ~ 130 cm LV 90 LM 110 LL 120 LE 130 XLV 90 XLS 100 XLS 100 XLM 110 XLL 120		ML 120
LARGE ~ 120 cm X-LARGE ~ 130 cm LL 120 LL 120 LE 130 XLV 90 XLS 100 XLS 100 XLM 110 XLL 120		ME 130
LARGE LM 110 ~ 120 cm LL 120 LE 130 XLV 90 X-LARGE XLS 100 ~ 130 cm XLL 120		LV 90
~ 120 cm LM 110 LL 120 LE 130 XLV 90 XLS 100 XLM 110 XLL 120		LS 100
LL 120 LE 130 XLV 90 XLS 100 XLM 110 XLL 120		LM 110
X-LARGE XLV 90 X-LARGE XLS 100 XLM 110 XLL 120		LL 120
X-LARGE XLS 100 ~ 130 cm XLM 110 XLL 120		LE 130
~ 130 cm XLM 110 XLL 120		XLV 90
~ 130 cm XLM 110 XLL 120	-	XLS 100
XLL 120		XLM 110
XLE 130		XLL 120
		XLE 130

This wraparound belongs to all those, that don't want to have a belt outside.

The interior belt is fixed in and cannot be removed.

431-LSW_COM



	Length
	XSV 90
X-SMALL	XSS 100
	XSM 110
< 99 cm	SXL 120
	XSE 130
	SV 90
CMALL	SS 100
SMALL ~ 100 cm	SM 110
	SL 120
	SE 130
	MV 90
MEDIUM	MS 100
$\sim 110 \text{ cm}$	MM 110
	ML 120
	ME 130
	LV 90
LARGE	LS 100
$\sim 120 \text{ cm}$	LM 110
~ 120 cm	LL 120
	LE 130
X-LARGE ~ 130 cm	XLV 90
	XLS 100
	XLM 110
	XLL 120
	XLE 130

Thanks to the multi-adjustable Stretch belt the new coat apron 431-LSW_COM provides optimal weight distribution on shoulder and hip. Comfortable and quick to put on and take off.

The model has the best comfort when standing, moving or sitting.

Options: Additional arm sleeves (left or right) for cardiologists (431-LSW_COM_KARD).

533-DSW





1		
	Length	
X-SMALL < 99 cm	XSV 90	
	XSS 100	
	XSM 110	
	SXL 120	
	XSE 130	
	SV 90	
014411	SS 100	
SMALL	SM 110	
~ 100 cm	SL 120	
	SE 130	
	MV 90	
	MS 100	
MEDIUM	MM 110	
~ 110 cm	ML 120	
	ME 130	
	LV 90	
	LS 100	
LARGE	LM 110	
~ 120 cm	LL 120	
	LE 130	
	XLV 90	
X-LARGE ~ 130 cm	XLS 100	
	XLM 110	
	XLL 120	
	XLE 130	

533-DSW coat aprons give wraparound protection from radiation during special examinations. They are characterised by the overlapping Velcro locking system on the shoulder. Thanks to its wide bottom edge the apron can be fastened to fit different sizes without any problems.

The height of the exterior belt is adjustable.

WIROMA by Lanz-Anliker AG	
	Clothing
Your specialist for radiation protection	Top and Skirt

437 & 537-Top / 538-Skirt





The overlapping Velcro lock at the front (437) and on the shoulder section (537) allow the top to be fastened comfortably and to fit well. Together with the 538-skirt it provides secure wraparound protection from radiation..

Optioin: Available in different lengths.

ТОР

Girth	Length
X-SMALL	Longar
	XSM 60
< 95 cm	V21/1 00
Umfang	
SMALL	
~ 100 cm	SM 60
Umfang	
MEDIUM	
~ 110 cm	MM 60
Umfang	
LARGE	
~ 120 cm	LM 60
Umfang	
X-LARGE	
~ 130 cm	XLM 60
Umfang	

SKIRT

Girth	Length
X-SMALL	
< 95 cm	XSM 60
Umfang	
SMALL	
~ 100 cm	SM 60
Umfang	
MEDIUM	
~ 110 cm	MM 60
Umfang	
LARGE	
~ 120 cm	LM 60
Umfang	
X-LARGE	
~ 130 cm	XLM 60
Umfang	

437 – TOP - KARD

!! Arm part also available with Wraparound apron 431-LSW_COM KARD **!!**

- Enlargement the shoulder
- Removable arm parts (left / right)
- Additional protection at longer work in x-ray beam



458-ARM KARD		
ltem-Nr.		Lead equivalent in X-Light material
45800200	One Size	0.25 mm Pb
45800300		0.35 mm Pb
45800400		0.50 mm Pb



The soft collar is used in areas where the user is exposed to X-ray radiation over a long period of time, like, for example: cardiac catheter, angiography, intramedullary bone nailing, femoral neck, etc.




- **4534 TS** Thyroid shield to protect thyroid glands. Its comfortable design guarantees protection and easy wearing. With Velcro strap / one size
- **4534 TSG** Can be fixed onto front or coat aprons or attached with press buttons.
- **4534 TSM** With magnetic fastening.

HYGIENIC MOLTON COVER



Composition: 68% molleton and 32 % PU (made to be antibacterial), More hygienic, washable at 40 °C



Gloves with opening for grip



One size for right and/or left hand. This article is available in different colours and lead values.

Available in: 0.35 / 0.50 / 0.75 mm Pb

Lead full-fingered gloves



Manufactured out of highly flexible lead vinyl. The exterior material is easy to maintain and made from artificial leather. The interior glove can be separated.

These gloves are available in one size, as a pair or individual items and in grey or blue: 0.50 mm Pb.

Surgical gloves



These sterilised disposable gloves have a rough surface.

Interventions suitable for radiology, orthopaedics, etc.

These gloves are available as a pair.



403 - X-ray head cover

Our Wiroma X-ray head cover is available in 3 different sizes, Small, Medium and Large including Velcro system for quick fastener. Available in lead-equivalent 0.35 mm and 0.50 mm Pb.

Quality: X-Light material.

The head cover can use in different departments:

- > Cardiology
- Diagnostic radiology
- > Angiography
- Radiology intervention
- ➢ etc.





Using radiation protection glasses considerably reduces the ocular lens' exposure to radiation. 90-99 % of scattered rays at 80 kV are already absorbed by a lead equivalent value of 0.50 mm. We offer different glasses designs to protect eyes.

These glasses can also be ordered with correcting glass, with single vision or bifocal lenses (additional charge). Please ask for our special form for the ophthalmologist to place these kinds of orders.











UNIVISTA TYP II 0.75mm Pb

Art. Nr. 31700901 Can be worn over glasses, the angle and length of temple can be adjusted. Weight: 96g

UNIUNITE 0.75mm Pb

White Art. Nr. 31700902 / Blue Art Nr. 31700903 Can be worn over glasses, good shielding with cheek/brwon protection. Available in wihte and blue. Weight: 76g

UNICRUISER 0.75mm Pb

Art. Nr. 31700125 Good shielding with this sporty plastic frame, correction in single vision and multi vision. Weightt: 74g

UNISURFER 0.75mm Pb

Art. Nr. 31700126 Good shielding withg this sporty plastic frame, high quality x-ray protective spectacles. Weight: 79g

UNICHARGER 0.75mm Pb

Art. Nr. 31700127 Softer nosepiece and temple ends, correction in simple vision and multi vision. Weight: 77g







MIRAGE / PERSOL 0.75 mm Pb Art. Nr. 463210 Correction in single vision, bifocal and multi vision. Weight: 92g



UNICLASS / URBAN 0.75 mm Pb Art. Nr. 462570 Good shielding with this sporty plastic frame, stlyisch and high quality. Weight: 72g





UNIPRO / CASUAL 0.75 mm Pb Art. Nr. 461830

High quality metal frame version easily adjusable, flexible and soft temple ends, correction in single visoon and multi vision. Weight: 71g

UNICAJA / CASUAL Light 0.75 mm Pb

Art. Nr. 461850 High quality metal frame version easily adjusable, correction in single vision and multi vision. Weight: 71g



417 – OPT



417-OPT for panoramicradiography	
Lead value	Girth x Length
0.25 mm Pb	60 x 60
0.35 mm Pb	60 x 60
0.50 mm Pb	60 x 60

Patient protection for panorama x-ray examination for children and adults.

420 – DTL

Our poncho is easy to put ouver the patient.

420-DTL for normal radiography	
Lead value	Girth x Length
0.35 mm Pb	60 x 88
0.50 mm Pb	60 x 88

418 – DTC



418-DTC for normal radiography	
Lead value	Girth x Length
0.25 mm Pb	60 x 88
0.35 mm Pb	60 x 88
0.50 mm Pb	60 x 88

Patients' apron with thyroid collar



419-DTL	
for normal radiography	
Lead value	Girth x Length
0.25 mm Pb	60 x 88
0.35 mm Pb	60 x 88
0.50 mm Pb	60 x 88

Patients' apron without thyroid collar



To protect patients' gonads during X-ray examinations for thorax and the upper abdominal area.



with plastic buckle

442-HSV with plastic buckle	
Umfang x Länge	
30 x 30	
40 x 30	
60 x 45	
60 x 60	





with Velcro strap



with hip wire

440-ASB with velcro strap	
Umfang x Länge	
0.50 mm Pb	30 x 30
	40 x 30
	60 x 45
	60 x 60

443-RSV With hip wire	
Umfang x Länge	
0.50 mm Pb	30 x 30
	40 x 30
	60 x 45
	60 x 60

Patient Protection

Childrens' aprons



488-Skirt:

For our smallest patients to protect the gonads during X-ray examinations of the thorax and upper abdominal area.

487-Top:

The overlapping Velcro strap makes this top easy to close and fit snugly on children.

483-AWC

Front apron for children. Cross-over Velcro strap with stretch over the hips.

480-FCW (LSW)

Coat apron for children made from highly flexible multi-layered protective material with very durable top fabric.

477-OPT

Children's dental apron. Protects patients during dental radiography, particularly on the dental chart.

All designs are available in lead equivalent values of 0.35 mm Pb or 0.50 mm Pb.



Motifs:

by Lanz-Anliker AG





All-round radiation protection with Velcro straps to cover and wrap patients. Shape made to order!



Lead vinyl mattress for covering patients during diagnosis and for use in therapy. Shape made to order.



Special wrap for patients 0.35 mm Pb / X-light or normal material Pad: 55 x 60 cm, Cover: 85 x 35 cm Adjustable for slight and larger patients.



Studies have revealed that the lens is apparently more sensitive to radiation than previously thought. According to estimations the point at which demonstrable opacity of the lens occurs is during only one short exposure at 0.5-2 Gy. The point at which visual impairment occurs (cataract) is at 5 Gy.

The CT Eye Shield has been specially developed to protect the lens.

- Clear reduction and optimisation of the dose during the CT scan and other X-ray procedures.
- The fitted spacer reduces artefacts to an absolute minimum.
- Hypoallergenic adhesive ensures that the shield consequently stays on.
- Material made from barium sulphate can be disposed of with normal waste without danger.

CT Thyroid Shield

The CT Thyroid Shield was designed to ensure a reduced dose on sensitive thyroid glands. It has a lead equivalent value of 0.125 mm Pb.




- A dose reduction of practically over 40 % during CT scans or general X-ray examinations.
- Over 60% reduction in the dose with transmissions of 50 kV to 150 kV.
- Material made from barium sulphate can be disposed of with normal waste without danger.



448 - MAMMO

For convering and protection of the left and the right breast during mammography.





5-Arm swivel bracket or 10-Arm swivel bracket



Out of stainless steel with moveable hangers (right and left). Lets you store your radiation protection wear while saving space and protecting the material. Includes fitting to make mounting it onto the wall easier.

With welded-on hooks.

Mobile suspension device





Made from sturdy stainless steel with moveable hangers. Lets you store your radiation protection wear while saving space, in a way that is both mobile and which protects the material.

With welded-on hooks.

for 5 aprons Dimensions: 60 x 55 x 165 cm (WxDxH)

for 10 aprons Dimensions: 88 x 70 x 155 cm (WxDxH) Made

Hanger

Wooden hanger Hanger made from solid beech wood. The defined shape of the hanger makes it suitable for all front and coat aprons, as well as tops and skirts.



Clothes hanger INOX

Hanger made from stainless steel. Fully welded. No moveable parts and struts. It is therefore very hygienic and easy to clean. Specially suited to operating theatres.

Protective curtain



Protective curtain for windows and room partitions, mobile design, made to measure. The curtains are available with a lead equivalent value of 0.25 mm - 1.00 mm Pb. We will gladly advise you on sizes and colours.

Possible uses:

- Gangways
- Entrance to control room
- Round curtain
- Shielding patients during X-ray procedures in the room
- Dividing different work places
- Covering up windows
- Double and multiple curtains

Mobile racks available in different sizes.